

ASTM D445, ASTM D446  
ASTM D2170, AASHTO T201

# BitUVisc viscometer

FULLY AUTOMATIC VISCOMETER FOR HIGHLY  
VISCOUS SAMPLES UP TO 120,000 mm<sup>2</sup>/s

"Finally an instrument that eliminates the need for the tedious, labour intensive manual loading, measurement and cleaning of high viscosity samples".



FULLY  
AUTOMATIC



HIGHLY  
ACCURATE



HANDLES  
EXTREME SAMPLES



EASY  
TO USE



DUPLO  
VISCOMETER TUBES

## FEATURES

- For high temperature, high viscosity applications
- Fully automatic measurement and cleaning
- Stand alone or PC control
- Preheated dual solvent cleaning

Omnitek introduces the newest member to the U-Visc family of fully automatic viscometer systems. Specifically developed to cover sample types which pose great challenges to automatic sampling, measurement and cleaning, because of their high pour point and viscosity.

### Flexible and reliable

The instrument was designed to be versatile and flexible, yet easy to use to cover the hardest applications. Different models are available with 1 or 2 viscometer tubes, allowing for the analysis of a wide range of viscosities and sample types.

Viscometer tubes are based on the well-known and proven Ubbelohde design. Using specially designed thermal sensors, these tubes allow the detection of fully opaque fluids without any

problems or limitations and eliminate the need to use the reverse flow tubes. The standard configuration offers viscosity measurement up to 150°C and sample preheating up to 200°C.

Sophisticated temperature control ensures that measurements are carried out well within the required temperature stability. The system features an innovative single position auto-sampling which provides integrated heating of the sample up to 200°C and uses preheated solvent to effectively clean the viscometer tube. The BitUVisc uses the special designed duplo viscometer tubes for unstable samples such as vacuum residues, additives, crude oils, waxes, heavy fuel oils, polymers, asphalts, etc., which typically show poor determinability when sampled repeatedly. These tubes have 2 measuring sections with an approximately similar tube constant. This allows for an actual duplo determination while only sampling once.

The instrument requires a PC for full operation, which is controlled through an advanced PC software application but it can also be controlled as a stand-alone unit through the responsive color touchscreen interface. PC software allows the operator to specify optimized and fully customizable test methods for each individual sample with the click of a mouse.

ASTM D445, ASTM D446, ASTM D2170, AASHTO T201

## FULLY AUTOMATIC VISCOMETRY SYSTEMS FOR HIGHLY VISCOUS SAMPLES

### SPECIFICATIONS

Measuring range	Up to 120,000 mm <sup>2</sup> /s, up to 150°C
Temperature range	15 - 150°C *
Temperature stability	15°C to 100°C ± 0.01°C *, Up to 150°C, ± 0.03°C
Sample volume	12 ml
Number of solvents	2
Viscometer type	Ubbelohde based
Sensor type	Thermal
PC Control	Multiple instruments controlled with 1 PC
Data export	USB

\* For temperatures around ambient, an external cooling circulator is required



### AVAILABLE MODELS

	BitUVisc 110	BitUVisc 120
Nr. of baths	1	1
Tubes	1	2
Dimensions	38 x 62 x 78 cm	
Weight	54 kg	

### THE BitUVisc IS SUPPLIED WITH

20-fold viscometer tubes
High temperature bath oil
Reference standards
Advanced PC software
All parts required for standard operation




**OMNITEK**

Coenecoop 715 Tel : +31 (0)182-302990  
2741 PW Waddinxveen Fax : +31 (0)182-302999  
The Netherlands info@omnitek.nl

[www.omnitek.nl](http://www.omnitek.nl)

Your authorized distributor is:



Disclaimer - All specifications, images and information provided in this brochure are subject to change.